

## CM2500 Chamber:

- 25 Tire Chamber
- Push Button Hydraulic Electric Door Open / Close
- Electric model 220 volt or 460 volt, 3 phase
- Silent Mercury Contacts (on electric version)
- Steam model fincoil indirect heat plus radial heat for fast warm-up
- Honeywell 4-Point Recorder, digital anytime setting, password protected
- Omron PLC PLC controls the entire cycle, so if all the set points are not correct, the cure cycle will not begin
- 10" Touch Screen Control, diagnostic view
- Chamber wall temperature display
- Thermocouple plug-in for independent check
- 2 Stainless Steel Manifolds (inflate & envelope control pressure)
- Oil Filled Air Gauges
- 10 HP fan motor
- 25 Vacuum Hoses, Long
- 5 Inflate Hoses, Long, 25 Inlets
- 1 Pie Pan Deflector
- 1 Lift Out Bridge (monorail)
- 2 Upright Legs (stands)

